

# Patient Information Leaflet

Mynarix | Please read carefully before use

## HbA1c Immunoassay Control

A bi-level quality control material for monitoring precision and accuracy of HbA1c immunoassay systems.

## What this product is used for

Assayed quality control to verify HbA1c immunoassay performance. Not a calibrator or for treatment decisions.

## How to use it

Reconstitute with the specified volume of deionized water, mix gently and equilibrate before use.

## Important things to know

- For in-vitro diagnostic use by trained personnel only.
- Handle per standard laboratory biosafety practice.
- Do not use beyond stated open-vial or expiry stability.

## Storing this product

Store lyophilized product at 2-8 deg C. Reconstituted control is stable per the value sheet when capped at 2-8 deg C.

## Questions or concerns

Talk to your doctor or pharmacist. To report a suspected side effect, contact Mynarix at [safety@mynarix.com](mailto:safety@mynarix.com) or +1 817 256 2824.

Mynarix, 2400 Skyline Drive, Suite 120, Mesquite Industrial District, Dallas, TX 75149 | +1 817 256 2824  
| [hello@mynarix.com](mailto:hello@mynarix.com)