

Patient Information Leaflet

Mynarix | Please read carefully before use

Lipid Panel Calibrator

A multi-analyte calibrator for standardising lipid panel assays across the clinical chemistry workflow.

What this product is used for

Calibration of lipid panel assays on compatible clinical chemistry analysers.

How to use it

Use per the analyser calibration protocol; verify with quality control material after calibration.

Important things to know

- For in-vitro diagnostic use by trained personnel only.
- Do not use as a control material.
- Discard if turbidity or contamination is observed.

Storing this product

Store at 2-8 deg C. Do not freeze. Return to refrigeration promptly after use.

Questions or concerns

Talk to your doctor or pharmacist. To report a suspected side effect, contact Mynarix at safety@mynarix.com or +1 817 256 2824.

Mynarix, 2400 Skyline Drive, Suite 120, Mesquite Industrial District, Dallas, TX 75149 | +1 817 256 2824
| hello@mynarix.com