

Mynarix

PRESCRIBING INFORMATION

Lipid Panel Calibrator

Diagnostics | Calibrator | Rx only

A multi-analyte calibrator for standardising lipid panel assays across the clinical chemistry workflow.

1. Composition

Analytes: Cholesterol, HDL, LDL, Triglycerides

Presentation: 5 x 3.0 mL

Traceability: Reference-method aligned

Storage: 2-8 deg C

2. Indications

Calibration of lipid panel assays on compatible clinical chemistry analysers.

3. Dosage and Administration

Use per the analyser calibration protocol; verify with quality control material after calibration.

4. Mechanism of Action

Provides assigned analyte values traceable to recognised reference measurement procedures.

5. Warnings and Precautions

- For in-vitro diagnostic use by trained personnel only.
- Do not use as a control material.
- Discard if turbidity or contamination is observed.

6. Storage

Store at 2-8 deg C. Do not freeze. Return to refrigeration promptly after use.

7. Manufacturer

Mynarix, 2400 Skyline Drive, Suite 120, Mesquite Industrial District, Dallas, TX 75149
| +1 817 256 2824 | hello@mynarix.com

For informational purposes only and not a substitute for professional medical advice. Approved indications and formulations vary by market.